

Product Datasheet

Human UBE2D1 protein (orb392656)

Description

Recombinant human UBE2D1 protein is

and the target gene en...

Reactivity

Human

Endotoxins

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

Tested

ELISA, MS, SDS-PAGE, WB

Applications
Target

UBE2D1

Preservatives

Lyophilized from a 0.2 μM filtered solution of 50mM HEPES, 150mM NaCl, 2mM DTT, 10% Glycerin pH 7.5.

Form/Appearance

Lyophilized from a 0.2 μM filtered solution of 50mM HEPES, 150mM NaCl, 2mM DTT, 10% Glycerin pH 7.5.

Storage

Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.

Note

For research use only

Application notes

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Dissolve the lyophilized protein in 1X PBS. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

Protein Sequence

MSPILGYWKI KGLVQPTRLL LEYLEEKYEE
HLYERDEGDK WRNKKFELGL EFPNLPYYID
GDVKLTQSMA IIRYADKHN MLGGCPKERA
EISMLEGAVL DIRYGVSRIA YSKDFETLKV
DFLSKLPEML KMFEDRLCHK TYLNGDHVTH
PDFMLYDALD VVLYMDPMCL DAFPKLVCFK
KRIEAIQID KYLKSSKYIA WPLQGWQATF
GGGDHPPKSD LVPRGSMALK RIQKELSDLQ
RDPPAHCSAG PVGDDLFHWQ ATIMGPPDSA
YQGGVFFLTV HFPTDYPFKP PKIAFTTKIY
HPNINSNGSI CLDILRSQWS PALTVSKVLL
SICSLLCDPN PDDPLVPDIA QIYKSDKEY
NRHAREWTQK YAM

Purity

Greater than 95% as determined by reducing SDS-PAGE.

Biorbyt Ltd.

7 Signet Court, Swann's Road, Cambridge, CB5 8LA, United Kingdom

Email: info@biorbyt.com | Phone: +44 (0) 1223 859-353 | Fax: +44 (0)1223 280 240

Biorbyt LLC.

68 TW Alexander Drive
Research Triangle Park
Durham, North Carolina
27709. United States

Email: info@biorbyt.com | Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558

Expiration Date

6 months from date of receipt.