

Product Datasheet

Human UBE2H protein (orb392660)

Description

Recombinant human UBE2H protein is

and the target gene enco...

Reactivity

Human

Endotoxins

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

Tested

ELISA, MS, SDS-PAGE, WB

Applications
Target

UBE2H

Preservatives

Liquid in a 0.2 μM filtered solution of 50mM HEPES, 150mM NaCl, 2mM DTT, 10% Glycerol, pH 7.5.

Form/Appearance

Liquid in a 0.2 μM filtered solution of 50mM HEPES, 150mM NaCl, 2mM DTT, 10% Glycerol, pH 7.5.

Storage

Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Dissolve the lyophilized protein in 1X PBS. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

Protein Sequence

MSPILGYWKI KGLVQPTRL LLEYLEEKYEE
HLYERDEGDK WRNKKFELGL EFPNLPYYID
GDVKLTQSMA IIRYADKHN MLGGCPKERA
EISMLEGAVL DIRYGVSRIA YSKDFETLKV
DFLSKLPEML KMFEDRLCHK TYLNGDHVTH
PDFMLYDALD VVLYMDPMCL DAFPKLVCFK
KRIEAIQID KYLKSSKYIA WPLQGWQATF
GGGDHPPKSD LVPRGSHMSS PSPGKRRMDT
DVVKLIESKH EVTILGGLNE FVVKFYGPQG
TPYEGGVWKV RVDLPDKYPF KSPSIGFMNK
IFHPNIDEAS GTVCLDVINQ TWTALYDLTN
IFESFLPQLL AYPNPIDPLN GDAAAMYLHR
PEEYKQKIKE YIQKYATEEA LKEQEEGTGD
SSSESSMSDF SEDEAQDMEL

Purity

Greater than 95% as determined by SEC-HPLC and reducing SDS-PAGE.