

Product Datasheet

Human UBE2J2 protein (orb392662)

Description

Recombinant human UBE2J2 protein is

and the target gene en...

Reactivity

Human

Endotoxins

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

Tested

ELISA, MS, SDS-PAGE, WB

Applications

Target

UBE2J2

Preservatives

Liquid in a 0.2 μM filtered solution of 20mM PB, 150mM NaCl, 1mM DTT, pH 7.2.

Form/Appearance

Liquid in a 0.2 μM filtered solution of 20mM PB, 150mM NaCl, 1mM DTT, pH 7.2.

Storage

Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Dissolve the lyophilized protein in 1X PBS. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

Protein Sequence

MSPILGYWKI KGLVQPTRL L LEYLEEKYEE
HLYERDEGDK WRNKKFELGL EFPNLPYYID
GDVKLTQSMA IIRYADKHN MLGGCPKERA
EISMLEGAVL DIRYGVSRIA YSKDFETLKV
DFLSKLPEML KMFEDRLCHK TYLNGDHVTH
PDFMLYDALD VVLYMDPMCL DAFPKLVCFK
KRIEAIQID KYLKSSKYIA WPLQGWQATF
GGGDHPPKSD LVPRGSPEFH MEKGPTLGS
ETSDFTKRQL AVQSLAFNLK DKVFCELFPE
VVEIKQKQK AQDELSSRPQ TLPLPDVVPD
GETHLVQNGI QLLNGHAPGA VPNLAGLQQA NR

Purity

Greater than 95% as determined by reducing SDS-PAGE.

Source

E. coli

Uniprot ID

[Q8N2K1](#)

Entrez

[118424](#)