

Product Datasheet

Human UBE2K protein (orb392664)

Description

Recombinant human UBE2K protein is

and the target gene enc...

Reactivity

Human

Endotoxins

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

Tested

Applications

ELISA, MS, SDS-PAGE, WB

Target

UBE2K

Preservatives

Liquid in a 0.2 μM filtered solution of 20mM PB, 150mM NaCl, pH 7.4.

Form/Appearance

Liquid in a 0.2 μM filtered solution of 20mM PB, 150mM NaCl, pH 7.4.

Storage

Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Dissolve the lyophilized protein in 1X PBS. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

Protein Sequence

MGHHHHHHGS DSEVNQEAKP EVKPEVKPET
HINLKVSDGS SEIFFKIKKT TPLRRLEAF
AKRQKGEMDS LRFLYDGIRI QADQAPEDLD
MEDNDIIEAH REQIGGMANI AVQRIKREFK
EVLKSEETSK NQIKVDLVDE NFTELRGEIA
GPPDTPYEGG RYQLEIKIPE TYPFNPPKVR
FITKIWHPI SSVTGAICLD ILKDQWAAAM
TLRTVLLSLQ ALLAAAEPPD PQDAVVANQY
KQNPEMFQQT ARLWAHVYAG APVSSPEYTK
KIENLCAMGF DRNAVIVALS SKSWDVETAT
ELLLSN

Purity

Greater than 95% as determined by reducing SDS-PAGE.

Source

E. coli

Uniprot ID

P61086

Entrez

3093