

## Product Datasheet

### Rat VEGFA protein (orb392698)

**Description**

Recombinant rat VEGFA protein is produced by

gene encodin...

**Reactivity**

Rat

**Endotoxins**

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

**Tested**

ELISA, MS, SDS-PAGE, WB

**Applications**
**Target**

VEGFA

**Preservatives**

Lyophilized from a 0.2 μM filtered solution of PBS, pH 7.4.

**Form/Appearance**

Lyophilized from a 0.2 μM filtered solution of PBS, pH 7.4.

**Storage**

Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.

**Note**

For research use only

**Application notes**

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Dissolve the lyophilized protein in 1X PBS. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

**Protein Sequence**

APTTEGEQKT HEVVKFMDVY QRSYCRPIET  
LVDIFQEYPD EIEYIFKPSC VPLMRCAGCC  
NDEALECVPT SESNVTMQIM RIKPHQSQHI  
GEMSFLQHSR CECRPKKDRT KPNHCEPCS  
ERRKHLFVQD PQTCKCCKN TDSRCKARQL  
ELNERTCRCD KPRR

**Purity**

Greater than 95% as determined by reducing SDS-PAGE.

**Source**

P. pastoris

**Uniprot ID**
**P16612**
**Entrez**
**83785**
**Expiration Date**

6 months from date of receipt.