

Product Datasheet

Human WDYHV1 protein (orb392716)

Description

Recombinant human WDYHV1 protein is

and the target gene en...

Reactivity

Human

Endotoxins

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

Tested

Applications

ELISA, MS, SDS-PAGE, WB

Target

WDYHV1

Preservatives

Liquid in a 0.2 μM filtered solution of PBS, 100mM GSH, 1% TritonX-100, 15% glycerol, pH 7.4.

Form/Appearance

Liquid in a 0.2 μM filtered solution of PBS, 100mM GSH, 1% TritonX-100, 15% glycerol, pH 7.4.

Storage

Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Dissolve the lyophilized protein in 1X PBS. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

Protein Sequence

MSPILGYWKI KGLVQPTRL L LEYLEEKYEE
HLYERDEGDK WRNKKFELGL EFPNLPYYID
GDVKLTQSMA IIRYADKHN MLGGCPKERA
EISMLEGAVL DIRYGVSRIA YSKDFETLKV
DFLSKLPEML KMFEDRLCHK TYLNGDHVTH
PDFMLYDALD VVLYMDPMCL DAFPKLVCFK
KRIEAIQID KYLKSSKYIA WPLQGQWQATF
GGGDHPPKSD LVPRGSMEGN GPAAVHYQPA
SPPRDACVYS SCYCEENVWK LCEYIKNHDQ
YPLEECYAVF ISNERKMIPI WKQQARPGDG
PVIWDYHVVL LHVSSGGQSF IYDLDTVLFP
PCLFDTYVED AIKSDDDIHP QFRRKFRVIC
ADSYLKNFAS DRSHMKDSSG NWREPPPPYP
CIETGDSKMN LNDFISMDPK VGWGAVYTLS
EFTHRFGSKN C

Purity

Greater than 95% as determined by reducing SDS-PAGE.

Biorbyt Ltd.

7 Signet Court, Swann's Road, Cambridge, CB5 8LA, United Kingdom

Email: info@biorbyt.com | Phone: +44 (0) 1223 859-353 | Fax: +44 (0)1223 280 240

Biorbyt LLC.

68 TW Alexander Drive
Research Triangle Park
Durham, North Carolina
27709, United States

Email: info@biorbyt.com | Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558

Expiration Date

6 months from date of receipt.