

Product Datasheet

Human ZBTB9 protein (orb392730)

Description

Recombinant human ZBTB9 protein is

and the target gene enc...

Reactivity

Human

Endotoxins

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

Tested

Applications

ELISA, MS, SDS-PAGE, WB

Target

ZBTB9

Preservatives

Lyophilized from a 0.2 μM filtered solution of PBS, pH 7.4.

Form/Appearance

Lyophilized from a 0.2 μM filtered solution of PBS, pH 7.4.

Storage

Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.

Note

For research use only

Application notes

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Dissolve the lyophilized protein in 1X PBS. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

Protein Sequence

METPTPLPPV PASPTCNAP RTIQIEFPQH
SSSLLLESLNR HRLEGKFCDV SLLVQGRELR
AHKAVLAAAS PYFHDKLLLG DAPRLTLPSV
IEADAFEGLL QLIYSGRLRL PLDALPAHLL
VASGLQMWQV VQCSEILRE LETSGGGISA
RGGNSYHALL STTSSTGGWC IRSPFQTPV
QSSASTESPA STESPVGGEG SELGEVLQIQ
VEEEEEEEED DDDDEDQGSAT LSQTPQPQRV
SGVFPRPHGP HPLPMTATPR KLPEGESAPL
ELPAPPALPP KIFYIKQEPF EPKEEISGSG
TQPGGAKEET KVFSGGDTEG NGELGFLLPS
GPGPTSGGGG PSWKPVDLHG NEILSGGGGP
GGAGQAVHGP VKLGGTPPAD GKRFGLCLGK
RFAVKPKRDR HIMLTFSLRP FGCGICNKRF
KLKHHLTEHM KTHAGALHAC PHCGRRFRVH
ACFLRHRDLC KGQGWATAHW TYKLEHHHHH H

Purity

Greater than 95% as determined by reducing

Biorbyt Ltd.

7 Signet Court, Swann's Road, Cambridge, CB5 8LA, United Kingdom

Email: info@biorbyt.com | Phone: +44 (0) 1223 859-353 | Fax: +44 (0)1223 280 240

Biorbyt LLC.

68 TW Alexander Drive
Research Triangle Park
Durham, North Carolina
27709, United States

Email: info@biorbyt.com | Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558