

Product Datasheet

Human ZMYND19 protein (orb392735)

Description

Recombinant human ZMYND19 protein is

and the target gene e...

Reactivity

Human

Endotoxins

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

Tested

Applications

ELISA, MS, SDS-PAGE, WB

Target

ZMYND19

Preservatives

Lyophilized from a 0.2 μM filtered solution of 20mM PB, 150mM NaCl, pH 7.2.

Form/Appearance

Lyophilized from a 0.2 μM filtered solution of 20mM PB, 150mM NaCl, pH 7.2.

Storage

Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.

Note

For research use only

Application notes

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Dissolve the lyophilized protein in 1X PBS. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

Protein Sequence

MGSSHHHHHH SSGLVPRGSH MTDKLGIVR
LGRVAGKTKY TLIDEQDIPL VESYSFEARM
EVDADGNGAK IFAYAFDKNR GRGSGRLLHE
LLWERHRGGV APGFQVVHLN AVTVDNRLDN
LQLVPWGWRP KAEETSSKQR EQSLYWLAIQ
QLPTDPIEEQ FVLNVTRY NANGDVVEEE
ENSCTYYECH YPPCTVIEKQ LREFNICGRC
QVARYCGSQC QKDWPAHKK HCRERKRPFQ
HELEPER

Purity

Greater than 95% as determined by reducing SDS-PAGE.

Source

E. coli

Uniprot ID

[Q96E35](#)

Entrez

[116225](#)