

Product Datasheet

GSK3B (Phospho-S9) Blocking Peptide (orb392769)

Description	The peptide is used to block Anti-GSK3 beta (pS9) antibody reactivity.
Reactivity	Human, Mouse, Porcine, Rat, Sheep, Zebrafish
Tested Applications	Blocking
Target	GSK3B
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry. Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	P18266, Q9WV60, P49841
Entrez	84027, 56637, 2932
Expiration Date	12 months from date of receipt.