

Product Datasheet

IGF1R (Phospho-Y1161) Blocking Peptide (orb392780)

Description	The peptide is used to block Anti-IGF1 Receptor (pY1161) antibody reactivity.
Reactivity	Human, Mouse, Rat, Zebrafish
Tested Applications	Blocking
Target	IGF1R
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry. Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	P24062, Q60751, P15208, P08069, P15127, P06213
Entrez	3480, 16337, 25718, 16001
Expiration Date	12 months from date of receipt.