

Product Datasheet

PRKD1 (Phospho-S738/S742) Blocking Peptide (orb392793)

Description	The peptide is used to block Anti-PRKD1/2/3 (pS738/S742) antibody reactivity.
Reactivity	Human, Mouse, Porcine, Primate, Rat, Sheep
Tested Applications	Blocking
Target	PRKD1; PRKD2; PRKD3
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q94806 , Q5XIS9 , Q9WTQ1 , Q8K1Y2 , Q62101 , Q15139 , Q8BZ03 , Q9BZL6
Entrez	23683 , 75292 , 85421 , 18760 , 292658 , 101540 , 25865 , 5587
Expiration Date	12 months from date of receipt.