

Product Datasheet

MAPK14 (Phospho-Y323) Blocking Peptide (orb392803)

Description	The peptide is used to block Anti-p38 (pY323) antibody reactivity.
Reactivity	Human, Mouse, Rat, Sheep
Tested Applications	Blocking
Target	MAPK14
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q16539 , P47811 , P70618
Entrez	1432 , 26416
Expiration Date	12 months from date of receipt.