

Product Datasheet

SIRT2 Blocking Peptide (orb392864)

Description	The peptide is used to block Anti-SIRT2 antibody reactivity.
Reactivity	Human, Mouse, Primate, Rat
Tested Applications	Blocking
Target	SIRT2
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q8VDQ8, Q8IXJ6, Q5RJQ4
Entrez	22933, 64383, 361532
Expiration Date	12 months from date of receipt.