

Product Datasheet

Synaptotagmin Blocking Peptide (orb392875)

Description	The peptide is used to block Anti-Synaptotagmin antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	SYT1; SYT2
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	P29101, P46096, P21707, P46097, P21579, Q8N9I0
Entrez	25716, 127833, 24805, 6857, 20980, 20979
Expiration Date	12 months from date of receipt.