

Product Datasheet

Adenosine A3 Receptor Blocking Peptide (orb392899)

Description	The peptide is used to block Anti-Adenosine A3 Receptor antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	ADORA3
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	P28647, P33765, Q61618
Entrez	140, 25370, 100911796
Expiration Date	12 months from date of receipt.