

Product Datasheet

ATP5H Blocking Peptide (orb392918)

Description	The peptide is used to block Anti-ATP5H antibody reactivity.
Reactivity	Human, Rat
Tested Applications	Blocking
Target	ATP5H
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	O75947 , P31399
Entrez	641434 , 10476
Expiration Date	12 months from date of receipt.