

Product Datasheet

Vasopressin peptide (orb392922)

| | |
|----------------------------|--|
| Description | The peptide is used to block Anti-Vasopressin V3 Receptor antibody reactivity. |
| Reactivity | Human |
| Tested Applications | Blocking |
| Target | AVPR1B |
| Form/Appearance | Lyophilized powder |
| Storage | Shipped at 4°C. Store at -20°C for one year. |
| Note | For research use only |
| Application notes | The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours. |
| Purity | > 85% |
| Source | Synthetic |
| Uniprot ID | P47901 |
| Entrez | 553 |
| Expiration Date | 12 months from date of receipt. |