

Product Datasheet

CTNNB1 (Phospho-T41/S45) Blocking Peptide (orb392963)

Description	The peptide is used to block Anti-Beta-catenin (pT41/S45) antibody reactivity.
Reactivity	Human, Mouse, Porcine, Rabbit, Rat, Sheep, Zebrafish
Tested Applications	Blocking
Target	CTNNB1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	P35222, Q02248, Q9WU82
Entrez	1499, 12387, 84353
Expiration Date	12 months from date of receipt.