

Product Datasheet

CTNND1 (Phospho-Y228) Blocking Peptide (orb392964)

Description	The peptide is used to block Anti-Catenin delta 1 (pY228) antibody reactivity.
Reactivity	Human, Mouse, Primate
Tested Applications	Blocking
Target	CTNND1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry. Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	O60716, P30999
Entrez	12388, 1500
Expiration Date	12 months from date of receipt.