

Product Datasheet

FER (Phospho-Y402) Blocking Peptide (orb392986)

Description	The peptide is used to block Anti-c-FER (pY402) antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	FER
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	P09760, P16591, P70451
Entrez	2241, 14158, 301737
Expiration Date	12 months from date of receipt.