

Product Datasheet

Histone H3 Blocking Peptide (orb393017)

Description	The peptide is used to block Anti-Histone H3 antibody reactivity.
Reactivity	Human, Mouse, Rat, Zebrafish
Tested Applications	Blocking
Target	HIST1H3A
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	P68431 , P84244 , P84243 , P84245 , P68433 , Q71DI3 , Q6LED0 , P84228
Entrez	15078 , 8358 , 8352 , 8351 , 8355 , 8357 , 8968 , 291159 , 100361558 , 8350 , 8353 , 8354 , 15077 , 8356 , 319152
Expiration Date	12 months from date of receipt.