

Product Datasheet

hnRNP K Blocking Peptide (orb393020)

Description	The peptide is used to block Anti-hnRNP K antibody reactivity.
Reactivity	Human, Mouse, Rabbit, Rat, Zebrafish
Tested Applications	Blocking
Target	HNRNPK
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	P61979, P61980, P61978
Entrez	117282, 3190, 15387
Expiration Date	12 months from date of receipt.