



Product Datasheet

Tubulin beta 3 antibody (orb393042)


Description	Rabbit polyclonal antibody to Tubulin beta 3.
Species/Host	Rabbit
Reactivity	A. thaliana, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IH, WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within the C-term region of human Beta-tubulin. The exact sequence is proprietary.
Target	TUBB3
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	Q4QRB4 , Q13509 , Q9ERD7
Entrez	246118 , 10381 , 22152
Dilution Range	WB: 1:500:1000, IHC-P: 1:100:200, IF/ICC: 1:100:500
Expiration Date	12 months from date of receipt.



Western blot analysis of Beta-tubulin ex...



Immunohistochemical analysis of Beta-tub...



Immunofluorescent analysis of Beta-tubul...