



Product Datasheet

GATA1 (Phospho-S142) antibody (orb393059)

Description	Rabbit polyclonal antibody to GATA1 (Phospho-S142)
Species/Host	Rabbit
Reactivity	Human, Mouse, Porcine, Primate, Rat
Conjugation	Unconjugated
Tested Applications	IH, IP, WB
Immunogen	KLH-conjugated synthetic phosphopeptide corresponding to residues surrounding S142 of human GATA1 protein. The exact sequence is proprietary.
Target	GATA1
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	P17679, P43429, P15976
Entrez	100911167, 25172, 2623, 14460
Dilution Range	WB: 1:500:1000, IHC-P: 1:100:200, IP: 1:10:100
Expiration Date	12 months from date of receipt.



Western blot analysis of GATA1 (Phospho-...



Immunohistochemical analysis of GATA1 (P...