


Product Datasheet

IGLC1 antibody (orb393079)

Description	Rabbit polyclonal antibody to IGLC1.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within the N-term region of human Ig lambda C. The exact sequence is proprietary.
Target	IGLC1; IGLC2; IGLC3; IGLC6; IGLC7
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	P0CG05 , P0CG06 , A0M8Q6 , P0CF74 , P0CG04 , P15814
Entrez	3543
Dilution Range	WB: 1:500:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of Ig lambda C exp...