

Product Datasheet

PRKD1 (Phospho-S738/S742) antibody (orb393091)

Description

Rabbit polyclonal antibody to PRKD1

Species/Host

Rabbit

Reactivity

Human, Mouse, Porcine, Primate, Rat, Sheep

Conjugation

Unconjugated

Tested

Applications

IH, WB

Immunogen

KLH-conjugated synthetic phosphopeptide corresponding to residues surrounding S738/S742 of human PRKD1/2/3 protein. The exact sequence is proprietary.

Target

PRKD1; PRKD2; PRKD3

Preservatives

Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.

Form/Appearance

Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.

Storage

Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.

Note

For research use only

Clonality

Polyclonal

Source

Rabbit

Uniprot ID

[Q9WTQ1](#), [Q8K1Y2](#), [Q15139](#), [Q5XIS9](#), [Q62101](#), [Q8BZ03](#), [Q9BZL6](#), [O94806](#)

Entrez


[85421](#), [25865](#), [18760](#), [292658](#), [75292](#), [101540](#), [5587](#), [23683](#)

Dilution Range


WB: 1:500:1000, IHC-P: 1:100:200

Expiration Date

12 months from date of receipt.



Western blot analysis of PRKD1/2/3 (Phos...)



Immunohistochemical analysis of PRKD1/2/...