


Product Datasheet

MAP3K1 antibody (orb393098)

Description	Rabbit polyclonal antibody to MAP3K1.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within the N-term region of human MEKK1. The exact sequence is proprietary.
Target	MAP3K1
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	Q13233
Entrez	4214
Dilution Range	WB: 1:500:1000, IF/ICC: 1:100:500
Expiration Date	12 months from date of receipt.



Western blot analysis of MEKK1 expression...