



Product Datasheet

MED1 (Phospho-T1457) antibody (orb393103)


Description	Rabbit polyclonal antibody to MED1 (Phospho-T1457)
Species/Host	Rabbit
Reactivity	Human, Mouse, Primate, Rat
Conjugation	Unconjugated
Tested Applications	ChIP, IF, IH, IP, WB
Immunogen	KLH-conjugated synthetic phosphopeptide corresponding to residues surrounding T1457 of human TRAP220 protein. The exact sequence is proprietary.
Target	MED1
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	Q925J9 , Q15648
Entrez	19014 , 5469
Dilution Range	WB: 1:500:1000, IHC-P: 1:100:200, IF/ICC: 1:100:500, IP: 1:10:100, ChIP: Use at an assay dependent concentration
Expiration Date	12 months from date of receipt.



Western blot analysis
of TRAP220
(Phospho-T1457)



Immunohistochemical
analysis of TRAP220



Immunofluorescent
analysis of TRAP220
(Phospho-T1457)