



Product Datasheet

PPP1R1B (Phospho-T34) antibody (orb393135)


Description	Rabbit polyclonal antibody to PPP1R1B (Phospho-T34)
Species/Host	Rabbit
Reactivity	Human, Mouse, Porcine
Conjugation	Unconjugated
Tested Applications	IF, IH, WB
Immunogen	KLH-conjugated synthetic phosphopeptide corresponding to residues surrounding T34 of human DARPP32 protein. The exact sequence is proprietary.
Target	PPP1R1B
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	Q9UD71 , Q60829
Entrez	19049 , 84152
Dilution Range	WB: 1:500:1000, IHC-P: 1:100:200, IF/ICC: 1:100:500
Expiration Date	12 months from date of receipt.



Western blot analysis of DARPP32 (Phospho-T34)



Immunohistochemical analysis of DARPP32



Immunofluorescent analysis of DARPP32 (Phospho-T34)