



Product Datasheet

RUNX1 (Phospho-S276) antibody (orb393156)

Description	Rabbit polyclonal antibody to RUNX1 (Phospho-S276)
Species/Host	Rabbit
Reactivity	Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	IH, WB
Immunogen	KLH-conjugated synthetic phosphopeptide corresponding to residues surrounding S276 of human RUNX1 protein. The exact sequence is proprietary.
Target	RUNX1
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	Q03347 , Q63046 , Q01196
Entrez	861 , 50662
Dilution Range	WB: 1:500:1000, IHC-P: 1:100:200
Expiration Date	12 months from date of receipt.



Western blot analysis of RUNX1 (Phospho-S276) antibody.



Immunohistochemical analysis of RUNX1 (Phospho-S276) antibody.