




Product Datasheet

KLRC1 antibody (orb39316)

Description	Peptide affinity purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC-P, WB
Immunogen	Synthetic Peptide
Target	KLRC1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	26314
Uniprot ID	P26715
NCBI	NP_998822.1 , NP_002250.1 , NP_015567.1 , NP_998823.1
Dilution Range	IF/ICC: 1:25, FACS: 1:25, WB: 1:1000, IHC-P: 1:50-100
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of MDA-MB435...



Immunohistochemical staining of paraffin...