

Product Datasheet

alpha Synuclein (Phospho-Y133) antibody (orb393167)

Description

Rabbit polyclonal antibody to alpha

Species/Host

Rabbit

Reactivity

Human, Mouse, Porcine, Primate, Rabbit, Rat, Sheep

Conjugation

Unconjugated

Tested Applications

IH, WB

Immunogen

KLH-conjugated synthetic phosphopeptide corresponding to residues surrounding Y133 of human Alpha-synuclein protein. The exact sequence is proprietary.

Target

SNCA

Preservatives

Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.

Form/Appearance

Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.

Storage

Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.

Note

For research use only

Clonality

Polyclonal

Source

Rabbit

Uniprot ID

[P37377](#), [O55042](#), [P37840](#)

Entrez


[29219](#), [6622](#), [20617](#)

Dilution Range


WB: 1:500:1000, IHC-P: 1:100:200

Expiration Date

12 months from date of receipt.



Western blot analysis of Alpha-synuclein...



Immunohistochemical analysis of Alpha-sy...