




Product Datasheet

KLRC2 antibody (orb39317)

Description	Rabbit polyclonal antibody to KLRC2
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IHC-P, WB
Immunogen	Synthetic Peptide
Target	KLRC2
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	26159
Uniprot ID	P26717
NCBI	NP_002251.2
Dilution Range	WB: 1:1000, FACS: 1:10-50, IHC-P: 1:50-100
Expiration Date	12 months from date of receipt.



Western blot analysis of in HL-60 cell l...



Western blot analysis of in K562 cell li...