



Product Datasheet

TLN1 antibody (orb393175)

Description	Rabbit polyclonal antibody to TLN1.
Species/Host	Rabbit
Reactivity	Human, Mouse, Porcine
Conjugation	Unconjugated
Tested Applications	IH, WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within the C-term region of human Talin 1. The exact sequence is proprietary.
Target	TLN1
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	P26039, Q9Y490
Entrez	7094, 21894
Dilution Range	WB: 1:500:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of Talin 1 express...



Immunohistochemical analysis of Talin 1 ...