



Product Datasheet

ACO1 (Phospho-S711) antibody (orb393196)

Description	Rabbit polyclonal antibody to ACO1 (Phospho-S711)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rabbit, Rat, Zebrafish
Conjugation	Unconjugated
Tested Applications	IH, WB
Immunogen	KLH-conjugated synthetic phosphopeptide corresponding to residues surrounding S711 of human Aconitase 1 protein. The exact sequence is proprietary.
Target	ACO1
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	Q63270 , P21399 , P28271
Entrez	50655 , 48 , 11428
Dilution Range	WB: 1:500:1000, IHC-P: 1:100:200
Expiration Date	12 months from date of receipt.



Western blot analysis of Aconitase 1 (Ph...



Immunohistochemical analysis of Aconitas...