


Product Datasheet

AKT1 (Phospho-T308/309) antibody (orb393202)

Description	Rabbit polyclonal antibody to AKT1 (Phospho-T308/309)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH-conjugated synthetic phosphopeptide corresponding to residues surrounding T308/309 of human AKT1/2 protein. The exact sequence is proprietary.
Target	AKT1; AKT2
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	P31751, Q60823, P31750, P47196, P31749, P47197
Entrez	11652, 11651, 208, 207, 25233, 24185
Dilution Range	WB: 1:500:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of AKT1/2 (Phospho-T308/309) antibody.