



Product Datasheet

ATF2 (Phospho-S62) antibody (orb393213)

Description	Rabbit polyclonal antibody to ATF2 (Phospho-S62)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IH, IP, WB
Immunogen	KLH-conjugated synthetic phosphopeptide corresponding to residues surrounding S62 of human ATF2 protein. The exact sequence is proprietary.
Target	ATF2
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	P15336, P16951, Q00969
Entrez	1386, 102641666, 81647, 11909
Dilution Range	WB: 1:500:1000, IHC-P: 1:100:200, IP: 1:10:100
Expiration Date	12 months from date of receipt.



Western blot analysis of ATF2 (Phospho-S...



Immunohistochemical analysis of ATF2 (Ph...