



Product Datasheet

BCR (Phospho-Y360) antibody (orb393223)

Description	Rabbit polyclonal antibody to BCR (Phospho-Y360)
Species/Host	Rabbit
Reactivity	Human, Mouse, Primate, Rat
Conjugation	Unconjugated
Tested Applications	IF, IH, WB
Immunogen	KLH-conjugated synthetic phosphopeptide corresponding to residues surrounding Y360 of human BCR protein. The exact sequence is proprietary.
Target	BCR
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	P11274, Q6PAJ1
Entrez	110279, 613
Dilution Range	WB: 1:500:1000, IHC-P: 1:100:200, IF/ICC: 1:100:500
Expiration Date	12 months from date of receipt.



Western blot analysis of BCR (Phospho-Y360) antibody.



Immunohistochemical analysis of BCR (Phospho-Y360) antibody.



Immunofluorescent analysis of BCR (Phospho-Y360) antibody.