



Product Datasheet

CHEK1 (Phospho-S345) antibody (orb393243)


Description	Rabbit polyclonal antibody to CHEK1 (Phospho-S345)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IH, WB
Immunogen	KLH-conjugated synthetic phosphopeptide corresponding to residues surrounding S345 of human CHK1 protein. The exact sequence is proprietary.
Target	CHEK1
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	Q91ZN7 , O14757 , O35280
Entrez	140583 , 1111 , 12649
Dilution Range	WB: 1:500:1000, IHC-P: 1:100:200, IF/ICC: 1:100:500
Expiration Date	12 months from date of receipt.



Western blot analysis of CHEK1 (Phospho-S...



Immunohistochemical analysis of CHEK1 (Ph...



Immunofluorescent analysis of CHEK1 (Phos...