


Product Datasheet

EGFR (Phospho-Y1092) antibody (orb393269)

Description	Rabbit polyclonal antibody to EGFR (Phospho-Y1092)
Species/Host	Rabbit
Reactivity	Human, Mouse, Primate, Rat
Conjugation	Unconjugated
Tested Applications	IH, WB
Immunogen	KLH-conjugated synthetic phosphopeptide corresponding to residues surrounding Y1092 of human EGFR protein. The exact sequence is proprietary.
Target	EGFR
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	P00533, Q01279
Entrez	13649, 1956
Dilution Range	WB: 1:500:1000, IHC-P: 1:100:200
Expiration Date	12 months from date of receipt.



Western blot analysis of EGFR (Phospho-Y1092) antibody.



Immunohistochemical analysis of EGFR (Phospho-Y1092) antibody.