



Product Datasheet

GFAP (Phospho-S38) antibody (orb393296)


Description	Rabbit polyclonal antibody to GFAP (Phospho-S38)
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IF, IH, WB
Immunogen	KLH-conjugated synthetic phosphopeptide corresponding to residues surrounding S38 of human GFAP protein. The exact sequence is proprietary.
Target	GFAP
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	P14136, P03995
Entrez	14580, 2670
Dilution Range	WB: 1:500:1000, IHC-P: 1:100:200, IF/ICC: 1:100:500
Expiration Date	12 months from date of receipt.



Western blot analysis of GFAP (Phospho-S38)



Immunohistochemical analysis of GFAP (Phospho-S38)



Immunofluorescent analysis of GFAP (Phospho-S38)