


Product Datasheet

Tubulin beta 3 antibody (orb393337)

Description	Mouse monoclonal antibody to Tubulin beta 3.
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence of human Beta3-tubulin. The exact sequence is proprietary.
Target	TUBB3
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Monoclonal
Source	Mouse
Uniprot ID	Q4QRB4 , Q13509 , Q9ERD7
Entrez	246118 , 10381 , 22152
Dilution Range	WB: 1:2000:5000
Expiration Date	12 months from date of receipt.



Western blot analysis of Beta3-tubulin e...