



Product Datasheet

Tubulin beta 3 antibody (orb393346)

Description	Rabbit polyclonal antibody to Tubulin beta 3.
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IH, WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence of human Beta3-tubulin. The exact sequence is proprietary.
Target	TUBB3
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	Q4QRB4 , Q13509 , Q9ERD7
Entrez	246118 , 10381 , 22152
Dilution Range	WB: 1:2000:5000, IHC-P: 1:200:500
Expiration Date	12 months from date of receipt.



Western blot analysis of Beta3-tubulin e...



Immunohistochemical analysis of Beta3-tu...