

Product Datasheet

CXCL11 peptide (orb393918)

Description	CXCL11 peptide
Conjugation	Unconjugated
Tested Applications	HPLC
Target	CXCL11
Form/Appearance	Lyophilized
Storage	Upon reconstitution, aliquot and store at -20°C or -80°C. for long term storage, Avoid repeated freeze/thaw cycles.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb13425 (CXCL11 antibody)
Uniprot ID	Q8R392, Q7TNL0
Expiration Date	6 months from date of receipt.