

Product Datasheet

VEGF165 peptide (orb393928)

Description	VEGF165 peptide
Conjugation	Unconjugated
Tested Applications	HPLC
Target	VEGF165
Form/Appearance	Lyophilized
Storage	Upon reconstitution, aliquot and store at -20°C or -80°C. for long term storage, Avoid repeated freeze/thaw cycles.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb181784 (VEGF165 antibody)
Expiration Date	6 months from date of receipt.