



Product Datasheet

KLC1 antibody (orb39409)


Description	KLC1 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Recombinant Human Kinesin light chain 1 protein (274-573AA)
Target	KLC1
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q07866
Dilution Range	WB:1:1000-5000, IHC-P:1:20-200
Expiration Date	12 months from date of receipt.



Western blot analysis of COLO 320(lane 1...



Immunohistochemical staining of human brain tissue showing brown staining.



Immunohistochemical staining of human placenta tissue showing brown staining.