



Product Datasheet

HNF1A antibody (orb39410)

Description	Purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Synthetic Peptide
Target	HNF1A
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	67356
Uniprot ID	P20823
NCBI	NP_000536.5
Dilution Range	WB: 1:1000, IF/ICC: 1:10-50
Expiration Date	12 months from date of receipt.



Confocal
immunofluorescent
analysis of H...