

Product Datasheet

Beta Actin antibody (orb394145)

Description	Rabbit polyclonal antibody to Beta Actin.
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat, Zebrafish
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC-P, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB: 1:1000-10000, IHC-P: 1:200, IF/ICC: 200
Expiration Date	12 months from date of receipt.