



Product Datasheet

acid phosphatase, prostate antibody (orb394662)

Description	Rabbit polyclonal antibody to acid phosphatase, prostate.
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	acid phosphatase, prostate
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	48 kDa
Uniprot ID	P15309
Dilution Range	WB:1:500-1:5000, IHC:1:20-1:200, IF/ICC:1:20-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of PC-3 cells using...



Immunohistochemical staining of paraffin...