



Product Datasheet

AFM antibody (orb394678)

Description	Mouse monoclonal antibody to AFM.
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, IP, WB
Immunogen	afamin
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	≥95% as determined by SDS-PAGE
MW	69 kDa, observed 85 kDa
Uniprot ID	P43652
Dilution Range	WB:1:200-1:2000, IP:1:200-1:2000, IHC:1:20-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of human blood tis...



Immunohistochemical staining of paraffin...